# remel

## RapID™ NF Plus System

### INTENDED USE

Remel RapID NF Plus System is a qualitative micromethod employing conventional and chromogenic substrates for the identification of medically important glucose-nonfermenting, gram-negative bacteria and other select glucose-fermenting, gram-negative bacteria not belonging to the family *Enterobacteriaceae*, which have been isolated from human clinical specimens. A complete listing of the organisms addressed by the RapID NF Plus System is provided in the RapID NF Plus Differential Chart.

### SUMMARY AND EXPLANATION

The RapID NF Plus System is comprised of (1) RapID NF Plus Panels and (2) RapID NF Plus Reagent. Each RapID NF Plus Panel has several reaction cavities molded into the periphery of a plastic disposable tray. Reaction cavities contain dehydrated reactants and the tray allows the simultaneous inoculation of each cavity with a predetermined amount of inoculum. A suspension of the test organism in RapID Inoculation Fluid is used as the inoculum which rehydrates and initiates test reactions. After incubation of the panel, each test cavity is examined for reactivity by noting the development of a color. In some cases, reagents must be added to the test cavities to provide a color change. The resulting pattern of positive and negative test scores is used as the basis for identification of the test isolate by comparison of results to reactivity patterns stored in an Electronic RapID Compendium (ERIC™) database or by use of the RapID NF Plus Differential Chart.

### PRINCIPLE

The tests used in the RapID NF Plus System are based upon the microbial degradation of specific substrates detected by various indicator systems. The reactions employed are a combination of conventional tests and single-substrate chromogenic tests and are described below in Table 1.

### REAGENTS\*

**RapID NF Plus Reagent** (provided with kit) (15 ml/Bottle)  
Reactive Ingredient per Liter:  
p-Dimethylaminocinnamaldehyde ..... 0.05 g

**RapID Inoculation Fluid** (R8325102, supplied separately) (1 ml/Tube)  
KCl ..... 6.0 g  
CaCl<sub>2</sub> ..... 0.5 g  
Demineralized Water ..... 1000.0 ml

**RapID Nitrate A Reagent** (R8309003, supplied separately) ..(15 ml/Bottle)  
Sulfanilic Acid ..... 8.0 g  
Glacial Acetic Acid ..... 280.0 ml  
Demineralized Water ..... 900.0 ml

**RapID Spot Indole Reagent** (R8309002, supplied separately)(15 ml/Bottle)  
p-Dimethylaminocinnamaldehyde ..... 10.0 g  
Hydrochloric Acid ..... 100.0 ml  
Demineralized Water ..... 900.0 ml

\*Adjusted as required to meet performance standards.

### PRECAUTIONS

This product is for *In Vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, media, and test panels after use. Directions should be read and followed carefully.

#### Caution!

1. RapID NF Plus Reagent is toxic and may cause harm to the environment. Harmful by inhalation, contact with skin or eyes, or if swallowed. May impair fertility or cause harm to unborn child.
2. RapID Nitrate A Reagent and RapID Spot Indole Reagent may cause irritation to skin, eyes, and respiratory system.
3. Refer to Material Safety Data Sheet for detailed information on reagent chemicals.

### STORAGE

RapID NF Plus System, RapID Nitrate A, and RapID Spot Indole Reagents should be stored in their original containers at 2-8°C until used. Allow products to equilibrate to room temperature before use. DO NOT interchange reagents among different RapID systems. Remove only the number of panels necessary for testing. Reseal the plastic pouch and promptly return to 2-8°C. Panels must be used the same day they are removed from storage. RapID Inoculation Fluid should be stored in its original container at room temperature (20-25°C) until used.

### PRODUCT DETERIORATION

This product should not be used if (1) the color of the reagent has changed, (2) the expiration date has passed, (3) the plastic tray is broken or the lid is compromised, or (4) there are other signs of deterioration.

### SPECIMEN COLLECTION, STORAGE, AND TRANSPORT

Specimens should be collected and handled following recommended guidelines.<sup>11,12</sup>

**Table 1. Principles and Components of the RapID NF Plus System**

Cavity #	Test Code	Reactive Ingredient	Quantity	Principle	Bibliography #
<b>Before Reagent Addition:</b>					
1	ADH	Arginine	1.0%	Hydrolysis of arginine releases basic products which raise the pH and change the indicator.	1-3
2	TRD	Aliphatic thiol	0.2%	Utilization of the substrate lowers the pH and changes the indicator.	3
3	EST	Triglyceride	1.0%	Hydrolysis of the lipid releases fatty acids which lower the pH and change the indicator.	1-4
4	PHS	p-Nitrophenyl-phosphoester	0.1%	Enzymatic hydrolysis of the colorless aryl-substituted glycoside or phosphoester releases yellow σ- or p-nitrophenol.	2,3,5,6
5	NAG	p-Nitrophenyl-N-acetyl-β,D-glucosaminide	0.1%		
6	αGLU	p-Nitrophenyl-α,D-glucoside	0.1%		
7	βGLU	p-Nitrophenyl-β,D-glucoside	0.1%		
8	ONPG	p-Nitrophenyl, β,D-galactoside	0.1%		
9	URE	Urea	0.25%	Hydrolysis of urea produces basic products which raise the pH and change the indicator.	1-3
10	GLU	Glucose	1.0%	Utilization of glucose lowers the pH and changes the indicator.	1-3
<b>After Reagent Addition:</b>					
4	PRO	Proline-β-naphthylamide	0.1%	Enzymatic hydrolysis of the arylamide substrate releases free β-naphthylamine which is detected with RapID NF Plus Reagent.	4, 6-10
5	PYR	Pyrrolidine-β-naphthylamide	0.1%		
6	GGT	γ-Glutamyl β-naphthylamide	0.1%		
7	TRY	Tryptophane β-naphthylamide	0.1%		
8	BANA	N-Benzyl-arginine-β-naphthylamide	0.1%		
9	IND	Tryptophane	0.4%	Utilization of tryptophane results in the formation of indole which is detected with RapID Spot Indole Reagent.	1-3
10	NO <sub>3</sub>	Sodium nitrate	1.0%	Utilization of nitrate ion results in the formation of nitrite which is detected with RapID Nitrate A Reagent.	1-3

### MATERIALS SUPPLIED

(1) 20 RapID NF Plus Panels, (2) 20 report forms, (3) RapID NF Plus Reagent (one plastic dropper-bottle containing reagent sufficient for 20 panels), (4) 2 chipboard incubation trays, (5) Instructions for use (IFU).

### MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Gram stain reagents, (7) Microscope slides, (8) Oxidase reagent, (9) Cotton swabs, (10) RapID Inoculation Fluid-1 ml (R8325102), (11) RapID Nitrate A Reagent (R8309003), (12) RapID Spot Indole Reagent (R8309002), (13) McFarland #1 and #3 turbidity standards or equivalents (R20411 & R20413), (14) Pipettes, (15) ERIC (Electronic RapID Compendium, R8323600).

### PROCEDURE

#### Inoculum Preparation:

1. Test organisms must be grown in pure culture and examined by Gram stain and oxidase prior to use in the system.

**Note:** The oxidase test should be interpreted with caution when using bacterial growth from differential agars that contain dyes which may interfere with interpretation.

2. Test organisms may be removed from a variety of selective and nonselective agar growth media. The following types of media are recommended: Tryptic Soy Agar (TSA) with or without 5% Sheep Blood; Nutrient Agar; Chocolate Agar; MacConkey Agar.

#### Notes:

- Some media containing or supplemented with mono- or disaccharides are not recommended since they may suppress glycolytic activity and reduce test selectivity.
- Plates used for preparing inoculum should preferably be 18-24 hours old. Slow-growing isolates may be tested using 48 hour plates.
- The use of media other than those recommended may compromise test performance.

3. Using a cotton swab or inoculating loop, suspend sufficient growth from the agar plate culture in RapID Inoculation Fluid (1 ml) to achieve a visual turbidity equal to at least a #1 McFarland turbidity standard or equivalent, but not in excess of a #3 McFarland standard or equivalent.


#### Notes:

- Suspensions significantly less turbid than a #1 McFarland standard will result in aberrant reactions.
- Bacterial suspensions that are more turbid than a #1 McFarland standard will not affect test performance and are recommended for stock cultures and quality control strains. However, suspensions prepared with a turbidity significantly greater than a #3 McFarland standard may compromise test performance.
- Suspensions should be mixed thoroughly and vortexed if required.
- Suspensions should be used within 15 minutes of preparation.


4. An agar plate may be inoculated for purity and any additional testing that may be required using a loopful of the test suspension from the inoculation fluid tube. Incubate the plate for 18-24 hours at 35-37°C.

#### Inoculation of RapID NF Plus Panels:

- Peel back the lid of the panel over the inoculation port by pulling the tab marked "Peel to Inoculate" up and to the left.
- Using a pipette, gently transfer the **entire** contents of the Inoculation Fluid tube into the upper right-hand corner of the panel. Reseal the inoculation port of the panel by pressing the peel-back tab back in place.
- After adding the test suspension, and while keeping the panel on a level surface, tilt the panel back away from the reaction cavities at approximately a 45-degree angle (see below).


- While tilted back, gently rock the panel from side to side to evenly distribute the inoculum along the rear baffles as illustrated below.

- While maintaining a level, horizontal position (best achieved by using the bench top against the reaction cavity bottoms), slowly tilt the panel forward toward the reaction cavities until the inoculum flows along the baffles into the reaction cavities (see below). This should evacuate all of the inoculum from the rear portion of the panel.

**Note:** If the panel is tilted too quickly, air may be trapped at the test cavity junction, restricting fluid movement.

- Return the panel to a level position. If necessary, gently tap the panel on the bench top to remove any air trapped in the cavities.

#### Notes:

- Examine the test cavities, which should appear bubble-free and uniformly filled. Slight irregularities in test cavity fills are acceptable and will not affect test performance. If the panel is grossly misfilled, a new panel should be inoculated and the misfilled panel discarded.
- Complete the inoculation of each panel receiving inoculation fluid before inoculating additional panels.
- Do not allow the inoculum to rest in the back portion of the panel for prolonged periods without completing the procedure.

#### Incubation of RapID NF Plus Panels:

Incubate inoculated panels at 35-37°C in a non-CO<sub>2</sub> incubator for 4 hours. For ease of handling, panels may be incubated in the chipboard incubation trays provided with the kit.

**Note:** If desired, after an incubation period of 4 hours, and prior to the addition of any reagents, RapID NF Plus Panels may be placed in the refrigerator (2-8°C) overnight for reading the next morning.

#### Scoring of RapID NF Plus Panels:

RapID NF Plus panels contain 10 reaction cavities that, in addition to oxidase, provide 18 test scores. Test cavities 4 through 10 are bifunctional, containing two separate tests in the same cavity. Bifunctional tests are first scored before the addition of reagent providing the first test result, and then the same cavity is scored again after the addition of reagent to provide the second test result. The bifunctional test cavities that require RapID NF Plus Reagent are indicated with the first test above the bar and the second test below the bar. Bifunctional test 9, which requires RapID Spot Indole Reagent, is designated by a line under the reagent-requiring test, IND. Bifunctional test 10, which requires RapID Nitrate A Reagent, is designated with a box drawn around the reagent-requiring test, NO<sub>3</sub>.

#### RapID NF Plus Panel Test Location

Cavity #	1	2	3	4	5	6	7	8	9	10
Test Code	ADH	TRD	EST	PHS	NAG	αGLU	βGLU	ONPG	URE	GLU
				PRO	PYR	GGT	TRY	BANA	<u>IND</u>	<u>NO<sub>3</sub></u>

- While firmly holding the RapID NF Plus panel on the benchtop, peel off the label lid over the reaction cavities by pulling the lower right hand tab up and to the left.
- Without the addition of any reagents, read and score cavities 1 (ADH) through 10 (GLU) from left to right using the interpretation guide presented in Table 2. Record the test scores in the appropriate boxes on the report form using the test code above the bar for bifunctional tests.

**Note:** Record the **color** of cavity 10 (GLU) in the space provided on the report form. Blue indicates alkalization, green indicates oxidation, and yellow indicates fermentation. This information may be useful as an additional characteristic in resolving probability overlaps.

- Add the following reagents to the cavities indicated:

## ENGLISH

- Add 2 drops of RapID NF Plus Reagent to cavities 4 (PRO) through 8 (BANA).
- Add 2 drops of RapID Spot Indole Reagent to cavity 9 (URE/IND).
- Add 2 drops of RapID Nitrate A Reagent to cavity 10 (GLU/NO<sub>3</sub>).

**Note:** Only RapID Spot Indole Reagent should be used. Kovacs' or Ehrlich's Indole reagent will not provide satisfactory results.

4. Allow at least 30 seconds but no more than 3 minutes for color development. Read and score cavities 4 through 10. Record the scores in the appropriate boxes of the report form using the test codes below the bar for bifunctional tests.
5. Record the oxidase reaction for the test isolate in the box provided on the report form.
6. Reference the microcode obtained on the report form in ERIC for the identification.

**Table 2. Interpretation of RapID NF Plus System Tests\***

Cavity #	Test Code	Reagent	Reaction		Comments
			Positive	Negative	
Before Reagent Addition:					
1	ADH	None	Red	Yellow or light orange	Only the development of a distinct red color should be scored as positive.
2	TRD	None	Yellow, gold, or yellow-orange	Red or orange	Development of yellow, gold, or yellow-orange throughout the cavity should be scored as positive. A yellowish layer may form on the top of the cavity. Gently shake the panel to mix and score as outlined above.
3	EST	None	Yellow, gold, yellow-orange, or light orange	Red or dark orange	Development of yellow, gold, yellow-orange, or light orange throughout the cavity should be scored as positive. A red layer may form on the top of the cavity. Gently shake the panel to mix and score as outlined above.
4	PHS	None	Yellow	Clear or light tan	Any development of a yellow color throughout the cavity should be scored as positive.
5	NAG				
6	αGLU				
7	βGLU				
8	ONPG				
9	URE	None	Red	Yellow, yellow-orange, or orange	Only the development of a distinct red color throughout the cavity should be scored as positive.
10	GLU	None	Yellow	Blue, blue-green, or green	Only a distinct yellow color throughout the cavity should be scored as positive. Record the cavity <b>color</b> in the space provided on the report form for reference if additional characteristics are needed for identification.
After Reagent Addition:					
4	PRO	RapID NF Plus Reagent	Purple, violet, red, dark orange, or dark pink	Clear, tan, light orange, or very pale pink	Only significant color development should be scored as positive. Pale shades of color should be scored as negative.
5	PYR				
6	GGT				
7	TRY				
8	BANA				
9	IND	RapID Spot Indole Reagent	Brown or black	Orange or red	Any development of a brown-black color should be scored as positive. Any other color should be scored as negative.
10	NO <sub>3</sub>	RapID Nitrate A Reagent	Red or orange	Clear, tan, or yellow	Any development of a red or orange color should be scored as positive. <b>Note:</b> RapID Nitrate B Reagent is NOT REQUIRED.

\*NOTE: Panels should be read by looking down through the reaction cavities against a white background.

## RESULTS AND RANGE OF EXPECTED VALUES

The RapID NF Plus Differential Chart illustrates the expected results for the RapID NF Plus System. Differential Chart results are expressed as a series of positive percentages for each system test. This information statistically supports the use of each test and provides the basis, through numerical coding of digital test results, for a probabilistic approach to the identification of the test isolate.

Identifications are made using individual test scores from RapID NF Plus panels in conjunction with other laboratory information (i.e., Gram stain, oxidase, growth on differential or selective media) to produce a pattern that statistically resembles known reactivity for taxa recorded in the RapID System database. These patterns are compared through the use of the RapID NF Plus Differential Chart, or by derivation of a microcode and the use of ERIC.

## QUALITY CONTROL

All lot numbers of the RapID NF Plus System have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality

control results are noted, patient results should not be reported. Table 3 lists expected results for the selected battery of test organisms.

## Notes:

- RapID reagent quality control is accomplished by obtaining the expected reactions for tests requiring the addition of the reagents (cavities 4-10).
- Organisms which have been repeatedly transferred on agar media for prolonged periods may provide aberrant results.
- Quality control strains should be stored frozen or lyophilized. Prior to use, quality control strains should be transferred 2-3 times from storage on an agar medium that is recommended for use with the RapID NF Plus System.
- Formulations, additives, and ingredients of culture media vary from manufacturer to manufacturer and may vary from batch to batch. As a result, culture media may influence constitutive enzymatic activity of designated quality control strains. If quality control strain results differ from the patterns indicated, a subculture onto medium from a different batch or from another manufacturer will often resolve quality control discrepancies.

**Table 3. Quality Control Chart for RapID NF Plus Panels**

Organism	ADH	TRD	EST	PHS	NAG	αGLU	βGLU	ONPG	URE	GLU	PRO	PYR	GGT	TRY	BANA	IND	NO <sub>3</sub>
<i>Acinetobacter baumannii</i> <sup>a</sup> ATCC® 19606	–	–	+	(–)	–	–	–	–	–	(+)	–	–	–	(–)	–	–	–
<i>Aeromonas hydrophila</i> <sup>b</sup> ATCC® 35654	+	+	+	+	+	–	+	+	(–)	+	+	(–)	(–)	–	(–)	+	+
<i>Elizabethkingia meningoseptica</i> <sup>c</sup> ATCC® 13253	(–)	(–)	V	(+)	+	+	(+)	(–)	–	–	+	+	+	+	+	+	–
<i>Oligella ureolytica</i> <sup>b</sup> ATCC® 43534	(–)	–	–	–	–	–	–	–	+	–	V	–	+	(–)	–	–	(+)

+, positive; –, negative; V, variable; (–), usually negative; (+), usually positive

<sup>a</sup>Previously designated as *Acinetobacter calcoaceticus*

<sup>b</sup>**Key indicator strains** demonstrate acceptable performance of the most labile substrate in the system and reactivity in a significant number of wells, according to Clinical and Laboratory Standards Institute recommendations for streamlined quality control.<sup>18</sup>

<sup>c</sup>Previously designated as *Flavobacterium meningosepticum*

## LIMITATIONS

1. The use of RapID NF Plus System and the interpretation of results requires the knowledge of a competent laboratorian who is trained in general microbiological methods and who judiciously makes use of training, experience, specimen information, and other pertinent procedures before reporting the identification obtained using this system.
2. Specimen source, oxidase reaction, Gram stain characteristics, and growth on selective agars should be considered when using the RapID NF Plus System.
3. The RapID NF Plus System must be used with pure cultures of test organisms. The use of mixed microbial populations or direct testing of clinical material without culture will result in aberrant results.
4. The RapID NF Plus System is designed for use with the taxa listed in the RapID NF Plus Differential Chart. The use of organisms not specifically listed may lead to misidentifications.
5. Expected values listed for RapID NF Plus System tests may differ from conventional test results or previously reported information.
6. The accuracy of the RapID NF Plus System is based upon the statistical use of a multiplicity of specially designed tests and an exclusive, proprietary database. The use of any single test found in the RapID NF Plus System to establish the identification of a test isolate is subject to the error inherent in that test alone.

## PERFORMANCE CHARACTERISTICS

The RapID NF Plus System performance characteristics have been established by laboratory testing of reference and stock cultures at Remel and by clinical evaluations using fresh clinical and stock isolates.<sup>13-17</sup>

## BIBLIOGRAPHY

1. Blazevic, D.J. and G.M. Ederer. 1975. Principles of Biochemical Tests in Diagnostic Microbiology. John Wiley & Sons, New York, NY.
2. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1999. Manual of Clinical Microbiology. 7<sup>th</sup> ed. ASM Press, Washington, D.C.
3. Holt J.G. and N.R. Krieg. 1984. Bergey's Manual of Systematic Bacteriology. Vol. 1. Williams and Wilkins, Baltimore, MD.
4. MacFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3<sup>rd</sup> ed. Lippincott Williams & Wilkins, Philadelphia, PA.
5. Gilardi, G.L., S. Hirsch, and M. Mandel. 1975. J. Clin. Microbiol. 1:384-389.
6. Humble, W.M., A. King, and I. Phillips. 1977. J. Clin. Pathol. 30:275-277.
7. Kilian, M. and P. Bulow. 1976. Acta. Pathol. Microbiol. Scand. B. 84:245-251.
8. Mulczyk, M. and A. Szewczuk. 1970. J. Gen. Microbiol. 61:9-13.
9. Nagatsu, T., M. Hino, H. Fuyamada, T. Hayakawa, S. Sakakibara, Y. Nakagowa, and T. Takemoto. 1976. Anal. Biochem. 74:466-476.
10. Westley, J.R., P.J. Anderson, V.A. Close, B. Halpern, and E.M. Lederberg. 1967. Appl. Microbiol. 15:822-825.
11. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry, and M.A. Pfaller. 2007. Manual of Clinical Microbiology. 9<sup>th</sup> ed. ASM Press, Washington, D.C.
12. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey and Scott's Diagnostic Microbiology. 12<sup>th</sup> ed. Mosby Elsevier, St. Louis, MO.
13. Eriqez, L.A., A.P. Jones, and N.E. Hodinka. 1991. Abstract C-217. Abstracts of the 91<sup>st</sup> General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
14. Kiska, D.L., A. Kerr, M.C. Jones, J.A. Caracciolo, B. Eskridge, M. Jordan, S. Miller, D. Hughes, N. King, and P.H. Gilligan. 1995. Abstract C-312. Abstracts of the 95<sup>th</sup> General Meeting of the American Society for Microbiology. ASM, Washington, D.C.

15. Kiska, D.L., A. Kerr, M.C. Jones, J.A. Caracciolo, B. Eskridge, M. Jordan, S. Miller, D. Hughes, N. King, and P.H. Gilligan. 1996. J. Clin. Microbiol. 34:886-891.
16. Kitch, T., M.R. Jacobs, and P.C. Appelbaum. 1991. Abstract C-215. Abstracts of the 91<sup>st</sup> General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
17. Kitch, T., M.R. Jacobs, and P.C. Appelbaum. 1992. J. Clin. Microbiol. 30:1267-1270.
18. Clinical and Laboratory Standards Institute (CLSI). 2008. Quality Control for Commercial Microbial Identification Systems; Approved Guideline. M50-A. CLSI, Wayne, PA.

## PACKAGING

REF R8311005, RapID NF Plus System..... 20 Tests/Kit

## Symbol Legend

<b>REF</b>	Catalog Number
<b>IVD</b>	In Vitro Diagnostic Medical Device
<b>LAB</b>	For Laboratory Use
	Consult Instructions for Use (IFU)
	Temperature Limitation (Storage Temp.)
<b>LOT</b>	Batch Code (Lot Number)
	Use By (Expiration Date)
<b>EC REP</b>	Authorized European Representative
	Manufacturer

RapID™ is a trademark of Thermo Fisher Scientific and its subsidiaries.  
ERIC™ is a trademark of Thermo Fisher Scientific and its subsidiaries.  
ATCC® is a registered trademark of American Type Culture Collection.

	12076 Santa Fe Drive Lenexa, KS 66215, USA <a href="http://www.remel.com">www.remel.com</a> , (800) 255-6730 International: (913) 888-0939
<b>EC REP</b>	Remel Europe Ltd. Clipper Boulevard West, Crossways Dartford, Kent, DA2 6PT, UK

For technical information contact your local distributor.

IFU 8311005, Revised September 4, 2013

Printed in U.S.A.

Rapid NF Plus Differential Chart

Organism	ADH	TRD	EST	PHS	NAG	αGLU	βGLU	ONPG	URE	GLU	PRO	PYR	GGT	TRY	BANA	IND	NO <sub>3</sub>	OXI
<i>Acinetobacter</i> spp.	2	2	96	13	3	0	0	0	7	50	5	4	1	18	7	0	0	0
<i>Actinobacillus ureae</i>	0	0	0	98	0	0	0	0	99	95	0	0	0	0	0	0	98	95
<i>Aeromonas caviae</i>	12	80	90	86	98	0	92	96	2	98	99	1	11	0	91	94	95	98
<i>Aeromonas hydrophila</i>	79	71	92	88	98	0	90	76	4	93	95	4	8	2	27	96	96	98
<i>Aeromonas veronii</i> - biogroup <i>sobria</i>	70	4	96	79	95	0	0	91	2	95	99	0	18	36	79	98	98	98
<i>Aeromonas veronii</i> - biogroup <i>veronii</i>	0	2	98	77	98	9	95	14	96	93	92	0	2	4	0	96	96	95
<i>Agrobacterium radiobacter</i>	26	6	1	77	98	96	96	26	94	1	92	98	31	76	92	0	88	98
<i>Alcaligenes faecalis</i>	0	56	3	79	0	2	0	0	2	0	4	0	2	11	9	0	2	99
<i>Alcaligenes piechaudii</i>	2	81	3	29	0	18	0	0	4	0	26	22	3	3	40	0	96	98
<i>Alcaligenes xylosoxidans</i>	6	96	4	86	0	2	0	0	9	0	57	98	4	6	1	0	98	98
<i>Bergeyella zoohelcum</i> <sup>a</sup>	87	22	2	76	0	0	0	0	95	0	78	0	93	98	84	87	0	98
<i>Bordetella bronchiseptica</i>	2	39	0	96	0	0	0	0	98	0	6	8	0	5	0	0	90	99
<i>Brevundimonas diminuta</i>	7	2	14	92	89	2	0	0	2	2	98	7	70	81	90	0	0	99
<i>Brevundimonas vesicularis</i>	0	8	40	71	5	92	88	0	0	0	90	0	9	90	98	0	7	99
<i>Burkholderia cepacia</i>	8	2	95	67	24	0	9	5	4	17	55	1	92	9	6	0	9	52
<i>Burkholderia gladioli</i>	0	0	99	67	81	0	0	0	0	0	90	90	57	87	0	0	2	2
<i>Burkholderia pseudomallei</i>	11	62	99	89	98	0	0	0	2	0	71	2	99	7	2	0	97	99
CDC IVC-2	1	11	86	33	0	4	0	0	99	0	12	15	2	2	0	0	0	98
CDC NO-1	0	0	88	2	0	0	0	0	0	0	0	0	0	2	0	0	98	0
<i>Chromobacterium violaceum</i>	96	4	30	79	82	0	0	0	0	79	4	20	99	0	11	6	56	90
<i>Elizabethkingia meningoseptica</i> <sup>b</sup>	2	12	28	88	98	96	91	24	3	2	90	96	95	95	96	96	2	98
<i>Comamonas acidovorans</i>	2	93	8	83	0	0	0	0	2	0	0	36	94	4	0	0	97	99
<i>Comamonas testosteroni</i>	0	73	72	62	0	0	0	0	2	0	0	0	39	0	0	0	93	98
<i>Flavobacterium breve</i>	0	40	0	87	0	0	0	0	7	0	97	98	96	95	97	93	0	98
<i>Flavobacterium IIb</i>	0	62	39	93	95	98	95	5	32	6	3	90	95	95	80	91	27	98
<i>Flavobacterium III</i>	79	0	4	30	98	96	97	69	20	0	0	92	95	96	95	92	0	98
<i>Kingella dentrificans</i>	0	0	4	0	0	0	0	0	0	93	98	0	0	0	0	0	92	98
<i>Kingella kingae</i>	0	0	0	88	0	0	0	0	0	81	11	0	0	69	1	0	0	99
<i>Methylobacterium</i> spp.	0	0	86	32	2	2	0	0	99	0	45	16	5	5	21	0	6	90
<i>Moraxella atlantae</i>	0	0	38	96	0	0	0	0	0	0	0	18	68	89	0	0	1	99
<i>Moraxella catarrhalis</i>	0	0	99	3	0	0	0	0	0	0	1	0	0	82	0	0	39	98
<i>Moraxella lacunata</i>	0	0	91	88	0	0	0	0	0	0	12	0	1	0	35	0	96	99
<i>Moraxella lincolni</i>	0	0	98	0	0	0	0	0	0	0	0	5	0	0	0	0	0	99
<i>Moraxella nonliquefaciens</i>	1	0	91	4	0	0	0	0	0	0	11	0	6	4	0	0	76	99
<i>Moraxella osloensis</i>	0	5	94	67	0	0	0	0	0	0	11	24	13	24	2	0	2	99
<i>Myroides odoratus</i> <sup>c</sup>	20	51	46	97	4	0	0	1	92	0	1	98	80	76	8	0	2	95
<i>Neisseria weaveri</i> / <i>elongata</i>	0	0	0	0	0	0	0	0	0	0	99	0	0	0	4	0	0	29
<i>Ochrobactrum anthropi</i> (Vd)	42	95	0	1	27	8	3	2	93	1	98	97	98	92	97	0	99	98
<i>Oligella ureolytica</i> (IVe)	4	8	1	2	0	0	0	0	98	0	71	4	99	12	0	0	16	99
<i>Oligella urethralis</i>	1	4	0	0	0	0	0	0	0	0	30	0	99	14	0	0	0	99
<i>Pasteurella aerogenes</i>	0	0	3	84	0	0	0	96	98	98	0	0	68	7	0	0	97	94
<i>Pasteurella haemolytica</i>	0	29	0	90	0	97	12	86	0	86	0	0	79	9	0	0	90	84
<i>Pasteurella multocida</i>	0	13	0	98	4	11	0	14	0	95	0	0	1	9	0	98	94	97
<i>Pasteurella pneumotropica</i>	0	0	0	98	0	82	0	83	94	90	0	0	90	0	0	90	95	92
<i>Plesiomonas shigelloides</i>	89	0	4	96	98	3	0	79	31	88	98	0	0	11	0	96	94	98
<i>Pseudomonas aeruginosa</i>	96	87	78	70	1	0	9	1	26	4	91	69	98	9	11	0	84	98
<i>Pseudomonas alcaligenes</i>	22	3	5	0	1	0	0	0	7	0	9	0	97	11	9	0	91	99
<i>Pseudomonas fluorescens</i> / <i>putida</i>	95	6	12	29	0	0	0	0	11	1	98	24	78	61	70	0	2	98
<i>Pseudomonas luteola</i> <sup>d</sup>	67	2	2	67	0	2	95	91	14	21	98	88	86	86	95	0	20	2
<i>Pseudomonas mendocina</i> (Vb-2)	98	2	24	84	0	0	0	0	5	0	99	0	98	14	69	0	97	99
<i>Pseudomonas oryzae</i> <sup>e</sup>	2	0	86	89	0	9	9	0	5	0	98	59	41	86	83	0	0	0
<i>Pseudomonas pseudoalcaligenes</i>	47	1	6	13	0	0	0	0	0	0	92	0	99	10	5	0	97	98
<i>Pseudomonas stutzeri</i> (Vb-1)	0	9	15	8	0	20	0	0	2	1	98	0	96	25	89	0	94	99
<i>Pseudomonas</i> Group 2	0	0	4	0	0	0	0	89	92	0	98	98	99	0	0	0	0	98
<i>Psychrobacter phenylpyruvicus</i> <sup>f</sup>	6	2	71	0	0	2	0	0	98	0	9	0	2	12	2	0	40	98
<i>Ralstonia pickettii</i> <sup>g</sup>	0	2	89	14	2	6	2	0	36	4	11	98	88	2	0	0	22	99
<i>Roseomonas</i> spp.	0	0	93	60	3	90	0	0	93	0	90	8	1	34	9	0	3	83
<i>Shewanella putrefaciens</i>	0	70	29	95	96	0	5	0	4	1	96	86	88	92	82	0	93	99
<i>Sphingobacterium multivorum</i> (Ilk-2)	4	5	7	18	97	92	96	88	95	3	1	90	73	90	14	0	9	89
<i>Sphingobacterium spiritivorum</i> (Ilk-3)	18	0	31	2	89	91	95	72	13	9	0	98	70	81	89	0	0	92
<i>Sphingomonas paucimobilis</i>	0	11	86	42	70	90	93	62	0	9	72	88	93	19	86	0	3	70
<i>Stenotrophomonas maltophilia</i>	18	86	88	91	37	24	91	2	30	4	98	5	97	38	72	0	14	4
<i>Suttonella indologenes</i>	13	15	49	90	0	0	0	0	0	95	7	0	92	40	28	99	0	98
<i>Vibrio alginolyticus</i>	0	0	87	95	95	2	3	0	0	90	98	48	26	22	97	89	95	98
<i>Vibrio cholerae</i>	0	2	28	98	98	0	2	88	0	98	9	0	5	16	0	98	98	98
<i>Vibrio damsela</i>	78	0	0	98	90	0	0	0	0	77	0	94	0	11	0	0	80	98
<i>Vibrio fluvialis</i> / <i>furnissii</i>	90	44	29	92	98	0	8	41	0	99	98	68	2	13	0	11	97	98
<i>Vibrio hollisae</i>	0	0	13	98	8	0	0	0	0	37	0	96	0	0	0	74	89	98
<i>Vibrio mimicus</i>	2	0	0	92	99	0	0	78	0	99	0	4	0	4	0	93	96	98
<i>Virbio parahaemolyticus</i>	2	11	88	96	84	0	0	2	9	98	98	56	0	18	98	97	98	97
<i>Vibrio vulnificus</i>	1	2	28	88	81	2	39	63	4	84	99	66	2	31	0	94	86	89
<i>Weeksella virosa</i> (Ilf)	0	98	5	34	0	0	0	2	0	0	0	98	5	98	96	78	0	98

<sup>a</sup>Previously designated as *Weeksella zoohelcum*<sup>b</sup>Previously designated as *Flavobacterium meningosepticum* and *Chryseobacterium meningosepticum*<sup>c</sup>Previously designated as *Flavobacterium odoratum*<sup>d</sup>Previously designated as *Chryseomonas luteola*<sup>e</sup>Previously designated as *Flavimonas oryzae*<sup>f</sup>Previously designated as *Moraxella phenylpyruvica*<sup>g</sup>Previously designated as *Burkholderia pickettii*

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- Eriqez, L.A., A.P. Jones, and N.E. Hodinka. 1991. Abstract C-217. Abstracts of the 91st General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Kiska, D.L., A. Kerr, M.C. Jones, J.A. Caracciolo, B. Eskridge, M. Jordan, S. Miller, D. Hughes, N. King, and P.H. Gilligan. 1995. Abstract C-312. Abstracts of the 95th General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Kiska, D.L., A. Kerr, M.C. Jones, J.A. Caracciolo, B. Eskridge, M. Jordan, S. Miller, D. Hughes, N. King, and P.H. Gilligan. 1996. J. Clin. Microbiol. 34:886-891.
- Kitch, T., M.R. Jacobs, and P.C. Appelbaum. 1991. Abstract C-215. Abstracts of the 91st General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Kitch, T., M.R. Jacobs, and P.C. Appelbaum. 1992. J. Clin. Microbiol. 30:1267-1270.
- Clinical and Laboratory Standards Institute (CLSI). 2008. Quality Control for Commercial Microbial Identification Systems; Approved Guideline. M50-A. CLSI, Wayne, PA.

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